

FRONT-LINE™

FIELD STERILIZER

When you need...

More **Capacity**

More **Simplicity**

More **Portability**

More **Reliability**



*Producing the world's only
steam sterilizer
designed for field use*



At FDI, our Vision is to "Accelerate healing with life-changing medical products and services."

The FRONT-LINE series of field sterilizers brings advanced sterilization capabilities closer to the point of injury.

- Surgical (Bookwalter® Retractor Sets)
- Dental (handpieces)
- Rigid containers / Sterilization trays

CAPACITY



Up to 25 lbs. of instruments

WATER



Less than 4 oz. of water per cycle

RUGGEDNESS



No transport case required

POWER




Dual power input
120/240 VAC., 1-PH, 50/60 Hz,
15 amps



Fort Defiance Industries, LLC
(865) 406-6348
sales@fortdefianceind.com

FDA CLEARED



FRONT-LINE™

FIELD STERILIZER

CAPACITY

- Accommodates a Bookwalter Retractor Set
- Processes up to 25 lbs. of surgical or dental instruments

SIMPLICITY

- Fully automated cycles
 - Immediate Use, Textiles, Wrapped Instruments, Dental Handpieces and Custom Cycles
 - 250°F, 270°F, and 275°F Cycle temperatures
- Easy to use touch screen with on-screen instructions
- Simple and reliable one-step door operation


PORTABILITY

- One integrated unit
- Rugged design with no transport case required
- Built-in handles

RELIABILITY

- Operating temperature: 35°F to 130°F
- Storage temperature: -65°F to 180°F
- Operating altitude: 0 to 10,000 ft
- Dynamic air removal technology

Model	Chamber	Maximum Load	Outside (LxWxH)	Weight	Capacity
FL120	12.0" dia. x 23.5" long	20 lbs.	31.5" x 20" x 22"	130 lbs.	Fits standard wrapped sterilization trays with outside dimensions of 10.5" x 20" x 3.5"
FL135	13.5" dia. x 24" long	25 lbs.	32" x 21.5" x 23.5"	140 lbs.	Fits full-length rigid sterilization containers with outside dimensions of 23.1" x 12.4" x 6.1"



Indications For Use Statement:

The Fort Defiance Industries FRONT-LINE Field Sterilizers are autoclaves designed for sterilizing heat- and moisture-stable medical, dental, and surgical materials, including wrapped and unwrapped, solid, porous, and hollow items (e.g., dental handpieces, suction pipes) in healthcare facilities.



Fort Defiance Industries, LLC
2411 Maremont Parkway | Loudon, TN 37774
(865) 406-6348 | sales@fortdefianceind.com

