

00805-61315
Transducer Removal Tool

00803-51150 - I.V.
Pole Spacer Tool

AMO Service Manual on CD
and literature for OPOMTWL,
OPOMTW2, and OPOR3019L

OPOR3019L - Laminar Flow
Phaco Tip 30 deg. (Qty=6),
with literature

S47927 - Diathermy
Service Tool Cable Assy

00805-90021
Diathermy
Calibration Load
Box Assy

00805-90020
Phaco Calibration
Load Box Assy

OPOMTWL
Laminar Flow
Metal Tip Wrench,
with literature

Abbott Medical Optics has
indicated that part number
OPOMTW2 - Metal Tip Wrench
Universal is now obsolete and non-
procureable. Mr. Lamar Reese
and Mr. Dennis Rumpel from
USAMMA Fort Detrick have
approved shipment of this kit
without this one item. If you have
any questions, please call either of
these individuals at (301) 619-4377
or (301) 619-1891, respectively.

S48981 - Vitrectomy Test
Plug

00803-90012 - I.V. Pole
Cable Pin Extractor Tool

00803-61351
Transducer
Calibration Tool

S263402 - Service
Phaco Housing
Straightening Tool

2103-5571 - Conn.
Coax Cable, BNC
male to Banana

00803-11801
Keyboard
Adapter Cable



SERVICE KIT, VITRECTOMY SYSTEM (ABBOTT MEDICAL OPTICS)

NSN: 6515-01-585-8857 P/N: AMOSK-001

FORT DEFIANCE INDUSTRIES CAGE CODE: 70784 (865) 408-0100

www.fortdefianceind.com