



# Reusable Bowie Dick Test Pack

## Save money and the environment:



The Steriscan Reusable Test Pack will generate significant cost and waste savings with each use.

## All-in-one solution:


The Steriscan Reusable Test Pack eliminates the need to carry multiple product lines and can be used as a Process Challenge Device or as an Air Removal Test.

## Easy to use and interpret:

The Steriscan Bowie Dick Indicator employs migrating front technology. No more guessing as to color change!




Unprocessed



Processed  
Reject



Processed  
Accept



U.S. Patent Nos. 7,740,802, 7,718,125 and 7,790,105. Other patents pending.