

**TOOL KIT MEDICAL
EQUIPMENT
MAINTENANCE AND
REPAIR BIOMEDICAL
UNIT**

**Support and Consumables
Handbook**

6545-01-662-0036

UA 8009

LIN W45197

REVISED

November 2017

U.S. Army Medical Materiel Agency

The information contained in this manual was prepared by the

U.S. Army Medical Materiel Agency
Office of PM, Medical Devices
Acquisition Information Division
693 Neiman Street
Fort Detrick, MD 21702-5001

Note: Unit Assemblage (UA) listings are different depending on when they were published. While the Line Item Number (LIN) for a particular set may remain the same from year to year, the National Stock Number (NSN) and the Unit Assemblage Code (UAC) of the set will change each time the UA is clinically reviewed and updated.

Different UA versions are provided on USAMMA website to assist you in viewing the current version, annotated with the Bill of Materials (BOM) Usage Code of 4. For your most accurate results, search for the UA using the NSN listed on your unit's property book listing. (Archived UAs will have a Usage Code of A.)

Support and Consumables Handbooks, as a rule, contain information for maintenance significant medical, dental, and veterinary equipment items and the associated accessories and consumables that are issued with the end item.

Click on the link below for a listing of everything included in a given UA. You need to enter the UAC that you are searching for.

https://app.usamma.amedd.army.mil/web_ua/ua_search.cfm

This manual is a living document and may contain inaccuracies due to changes in mission requirements or equipment updates.

To report corrections or additions please contact

Customer Relationship Management

by clicking above link or by calling
(301) 619-1288, DSN 343-1288

For additional equipment item information at the
USAMMA Medical Equipment Literature section
[click here.](#)

Acronyms used on Support and Consumables Handbook covers:

CSH	Combat Support Hospital
DES	Dental Equipment Set
DISS	Dental Instrumentation and Supply Set
DMS	Dental Materiel Set
FST	Forward Surgical Team
GYN	Gynecology
ICU	Intensive Care Unit
MES	Medical Equipment Set
MMS	Medical Materiel Set
ROC	Role of Care
TMDE	Test, Measurement, and Diagnostic Equipment
VES	Veterinary Equipment Set

**NOTE: CLICK ON TITLE
FOR DIRECT LINK**

TABLE OF CONTENTS

SECTION	DEVICE	NSN	PAGE
1	CASE TOOL AND ACCESSORIES W FOAM CUTOUTS	(J) 5140-01-482-2325	1

Please provide all comments on a
[Project Manager Medical Devices Comment Card.](#)

CASE TOOL AND ACCESSORIES W FOAM CUTOUTS

TOC

MANUFACTURER: HARDIGG CASES

Manufacturer Model# AL3627-1505-16045-200

Manufacturer's PHONE: (413) 665-2163 FAX: (413) 665-8330

Manufacturer's Website: <http://www.MilitaryCases.com>

AAC 'J' ITEM NSN 5140-01-482-2325

UAs: 8009

USAMMA POC: Customer Relations Management (301) 619-1288, DSN 343-1288

For current price information, refer to [MEDSILS](#) or [WebFLIS](#) (WebFLIS account required).

System Description: Tool case with foam cut-out recesses. Tools included in UA shown in Table.

NOTE: Tools not included with case.



CASE TOOL AND ACCESSORIES W FOAM CUTOUTS
PART # AL3627-1505-16045-200
NSN 5140014822325

**CASE TOOL AND ACCESSORIES W FOAM CUTOUTS
PHYSICAL & ENVIRONMENTAL SPECIFICATIONS**

Package Exterior Shipping Dimensions	Length (Inches)	0
	Width (Inches)	0
	Height (Inches)	0
	Cube (Cubic Feet)	0
	Weight	
Item Dimensions	Length (Inches)	38.87
	Width (Inches)	30
	Height (Inches)	23.67
	Cube (Cubic Feet)	15.9731
	Weight	50.30 lbs
Storage Temp		-20° to 140° F
Operating Temp		-20° to 140° F
Storage Humidity		
Operating Humidity		
Amps		
Watts		
Power Requirements	N/A	
Notes		

**CASE TOOL AND ACCESSORIES W FOAM CUTOUTS
ACCESSORIES & CONSUMABLES - STARTUP**

NSN	Notes	Nomenclature	MFR PN Hardigg Cases	DIST PN	QTY	UI	Unit Price	Total Price	Shelf Life Mon.
5140014822325	Tools found in UA listed in table below	Case Tool and Accessories w Foam Cutouts	AL3627- 1505- 16045-200		1	EA			N/A


**CASE TOOL AND ACCESSORIES W FOAM CUTOUTS
ACCESSORIES & CONSUMABLES**

NO RESUPPLY ITEMS

CASE TOOL AND ACCESSORIES W FOAM CUTOUTS**ITEM TABLES****Tools and Equipment in UA for NSN 5140-01-482-2325**

NSN	Nomenclature	PN	QTY	UI	MFG
5120-01-428-5657	Bit Set Screwdriver	J4739P	1	SE	Stanley Proto Industrial Tools
7920-01-479-0198	Brush Wire Scratch	44064	1	EA	Weiler Corp.
5210-01-477-3834	Caliper Vernier English & Metric 8"/200mm Range SS W/Case	125MEA-8/200	1	EA	L S Starrett Co.
7910-01-477-3836	Cleaner Vacuum Electric 3' Hose 18' Cord Nozzle & Crevice Tool	4Z906	1	EA	W W Grainger Inc
5130-01-459-6060	Drill Driver Elect Cordless 18 V Incl Charger 2 Batteries Case	DW959K-2	1	EA	Black and Decker
5133-01-478-7706	Drill Set Twist Component Quantity 29	B94.11M	1	SE	American Society of Mechanical
5120-01-335-2269	Extractor Set Screw Sz 1-6 Sq Shanks Plastic Pouch	52445	1	SE	Irwin Industrial Tool Co.
5110-00-247-2565	File Set Hand Qty 12 5 1/2in No 2 Cut Swiss Pattern	37398	1	EA	Cooper Industries Inc
4240-01-477-5380	Goggles Industrial D6 One Pcs Single Aperture Clear Anti-Fog	1VT68	1	EA	W W Grainger Inc
5120-01-428-5093	Hammer Hand 22 Carpenter's Nailing Curved-Claw 16 Ounce	51-505	1	EA	Stanley Proto Industrial Tools
6230-01-645-5157	Headlamp Tactical	AD-NOD-MED-BK	1	EA	Atlantic Diving Supply Inc
5120-00-293-3485	Pliers 10 In A16 Fence Single Head	B107.500	1	EA	American Society of Mechanical
5120-01-478-6295	Punch Set Center Solid Component Quantity 5 W/Plastic Case	S117PC	1	SE	L S Starrett Co.
5120-01-430-1734	Punch Set Drift Industrial Quality (Iq) Program Item	SPDL5K	1	SE	NMTC Inc
5120-01-430-5070	Riveter Blind Hand	200F	1	EA	NMTC Inc
5130-01-448-6702	Sander Disk Electric Portable 5" Hd Ac Only Random Orbit	DW420	1	EA	Black and Decker
5110-01-428-5422	Saw Hand Metal Cutting Blade Lg 10"-12" D-Hdl Grip	15-408	1	EA	Stanley Proto Industrial Tools
5120-00-851-4247	Screwdriver Bit Set and Handle		1	SE	
5120-01-112-9532	Screwdriver, Torque	B107.300	1	EA	American Society of Mechanical
5210-00-178-4227	Tape Measuring Stl 100' .375" W Orange Case Hand Crank	100L	1	EA	Cooper Industries Inc
6685-01-292-7873	Thermometer Self-Indicating Digitl Must Comply W/Prod Std R00044	FLUKE-51-2	1	EA	Fluke Corporation
5120-01-477-4762	Vise Machinist's V-Grooved Jaws W/Opening 2.625" 3.0" Jaw Widt	3ZH06	1	EA	General Tools Manufacturing
5120-01-398-7594	Wrench Adjustable Steel Open End Single Head Rack 12in Lg	944605000	1	EA	Sears Roebuck and Co
5120-01-336-5476	Wrench Pipe Aluminum Handle 457 Str Jaw 18in 2.5in Max Sgl End	31100	1	EA	Ridge Tool Co
5120-00-494-1911	Wrench Plier Curved Jaw With Cutter 8 1/2 In Lg	10WR	1	EA	Irwin Industrial Tool Co.
5120-01-482-7617	Wrench Set Combination 15 Piece Box & Open Ends 12 Pt Satin Finish	1200FASD	1	SE	Stanley Proto Industrial Tools
5120-01-430-4499	Wrench Set Socket 3/8" Drive 20 Pcs 6 Mm-19 Mm 12 Pt W/Hdl Ext	J52224	1	SE	Stanley Proto Industrial Tools
5120-01-478-5589	Wrench Set Socket 33 Components W/Metal Box	52138	1	SE	Stanley Proto Industrial Tools

APPENDIX
Additional Images for UA8009



CALIPER VERNIER ENGLISH & METERIC 8IN 200MM RANGE SS W CASE
PART # 125MEA-8/200
NSN 5210014773834

The 125 Vernier Caliper includes metric/inch readings, 0-6"/150mm caliper in plastic case. 200mm/8" Range, .02mm/.001" Grad.



**CLEANER VACUUM ELECTRIC 3FT HOSE 18FT CORD NOZZLE & CREVICE TOOL
PART # 4Z906
NSN 7910014773836**

This lightweight, dry pickup-only vacuum features a single-stage motor for powerful pickup and a dent-resistant, rustproof tank to withstand tough use



**DRILL DRIVER ELECT CORDLESS 18 V INCL CHARGER 2 BATTERIES CASE
PART # DW995K-2
NSN 5130014596060**



DRILL SET TWIST COMPONENT QUANTITY 29
PART # 99802
NSN 5133014787706

Contains 1 ea fractional sizes 1/16 to 1/2 in. In 1/64 in. Increments



FILE SET HAND QTY 12 5 .50IN NO 2 CUT SWISS PATTERN
PART # 37398
NSN 5110002472565

1 ea barrette, equaling, flat, marking, joint, knife, half round, round, slitting, square, three square and crossing needle files; 5-1/2 in. Lg; no. 2 cut



SANDER DISK ELECTRIC PORTABLE 5" HD AC ONLY RANDOM ORBIT
PART # DW420
NSN 5130014486702

Palm grip; random orbital; 2 amp; 12000 opm; pressure sensitive adhesive pad; dust sealed switch; double insulated



SCREWDRIVER, TORQUE
PART # TQSC6A
NSN 5120011129532



**THERMOMETER SELF-INDICATING DIGITAL MUST COMPLY W PROD STD R00044
PART # FLUKE-51-2
NSN 6685012927873**

Item must comply with requirements of defense supply center columbus production standard no. Lo1533



**WISE MACHINIST'S V-GROOVED JAWS W OPENING 2.625IN 3.0IN JAW WIDT
PART # 3ZH06
NSN 5120014774762**

Swivel vacuum base; v-grooved jaws; folds for storage; triangular base; 6-1/4 in. Base height; also furnished with slide-on rubber jaws for fragile work



**WRENCH SET COMBINATION 15 PIECE BOX & OPEN ENDS 12 PT SATIN FINISH
PART # J1200FASD
NSN 5120014827617**

12 point; sizes 5/16, 3/8, 7/16, 1/2, 9/16, 5/8, 11/16, 3/4, 13/16, 7/8, 15/16, 1, 1-1/16, 1-1/8 and 1-1/4 in.



Contains 1 each chrome plated, 12 point, 3/8 inch female square drive, standard length sockets, metric sizes 6, 7, 8, 9, 10, 11, 12, 13, 14, 15, 16, 17, 18, and 19 mm; 3 and 12 inch extensions; flex handle, speeder handle, ratchet handle, and universal joint.

**WRENCH SET SOCKET 0.375IN DRIVE 20 PCS 6 MM-19 MM 12 PT W HDL EXT
PART # J52224
NSN 5120014304499**



Contains 1 each chrome plated, 12 point, 3/8 inch female square drive, standard length, fractional sizes 3/8, 7/16, 1/2, 9/16, 5/8, 11/16, 3/4, 13/16, and 7/8 inch; deepwell sizes 3/8, 7/16, 1/2, 9/16, 5/8, 11/16, 3/4, 13/16, and 7/8 inch; universal sizes 7/16, 1/2, 9/16, 5/8, and 3/4 inch; 8 point, standard length sizes 1/4, 5/16, and 3/8 inch; includes ratchet, 3, 6, and 12 inch extensions; flex handle; speeder handle and universal joint.

**WRENCH SET SOCKET 33 COMPONENTS W METAL BOX
PART # 52138
NSN 5120014785589**